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TEAD Inhibitor Screening Platform (VitroX TEAD)

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
Selectivity and Efficacy Profile of TEAD Inhibitors

Biosensors rapidly evaluate drug efficacy against TEAD isoforms (TEAD1-4) within 5 minutes using microplate reader analysis.

- (A) SDX-301 exhibits broad inhibitory effects across YAP-TEAD interactions with potent activity specifically against TEAD1 and TEAD4.
 (B) SDX-301 exhibits broad inhibitory effects across TAZ-TEAD interactions with potent activity specifically against TEAD1.
 (C) VitroX TEAD provides quantitative IC₅₀ and R² values for all YAP/TAZ-TEAD isoform interactions.

Live-cell TEAD Inhibitor Screening Platform (LiveX TEAD)

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
Selectivity and efficiency profile of TEAD inhibitor across TEAD isoforms in live cells

The efficacy of TEAD inhibitors was evaluated in real-time in HEK293T cells expressing TEAD biosensors. After seeding the cells, drugs were diluted at various concentrations and applied, with measurements performed at desired time points using a microplate reader.

(A,B) The interaction between YAP/TAZ and TEAD isoforms was assessed at various concentrations and normalized against DMSO. IAG933 exhibited potent inhibitory activity against TEAD1 and TEAD4, but not against the others.

(C) Utilizing the LiveX TEAD platform enables the differentiation of drug specificity and assessment of drug persistence through real-time analysis.