WizLog[®]

Electronic Laboratory Notebook

2026-08-31

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

Composition

Ingredient	Compound	Weight (g)
PVC Resin	PV4310 () Suspension Grade (X) Emulsion Grade () Bulk Polymerised () Copolymer	5,000
Plasticiser	DBP/DOP/DIOP/DNP/DOTP	1,200
Stabiliser	Baerlocher Bearopan/Advapak NEO/Dibasic Lead Phosphite/Dibasic Lead Stearate/Tribasic Lead Stearate	200
Lubricant	Calcium Stearate/Stearic Acid	
Filler	Ground CaCO ₃ /Precipitated CaCO ₃	1,500
Pigment	DuPont/Tioxide	

Processing


1. Weigh all the ingredients and mix them except for PVC resin to a masterbatch.
2. Feed the PVC resin and the masterbatch to a compounding extruder and blend them to a compound with conditions below:
 - Temperature = 120°C
 - Screw Speed = 300 rpm
3. Cool the blended compound to 60°C and pelletize it to a suitable size.

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- Experimental Goal: Compare quantification of the 52 carbamates in pure water vs. in waste water
- List of Analyte:


Carbamates Analogs	Thio-Carbamates	Urea Analogs	Complexes & Chelators
Aldicarb Sulfoxide	Butylate	Bramacil	Febam (Fe+3 DMDTC)
Aldicarb Sulfone	Diallate	Chlororoxuran	Metam Na (MMDTC anion)
Aldicarb	DPTC	Diuran	DMDTC (Anion)
Carbaryl	Molinate	Fenuron	OH MMDTC (anion)
Carbofuran	Tillam	Monuron	Ziram (Zn-2 DMDTC)
3OH-Carbofuran	Triallate	Neburon	
Methiocarb	Vernolate	Siduron	
Methomyl		Tebuthiuron	
1-Naphthol			
- System Set-up: Alliance HT 2795/ZQ/PDA990
- Software: QuanOptimize in MassLynx 4.0
- Experimental Conditions:
 - HPLC: Symmetry C8 2.1x150mm, 35µm, 0.3mL/min, 40°C
 - Mobile Phase:
 - A: 10mM ammonium acetate in water at pH 5.0
 - B: 10mM ammonium acetate in acetonitrile
 - Ionization: ESI+
 - Capillary Voltage: 3.5kV
 - Source Temperature: 140°C
 - Desolvation Temperature: 350°C
- Data & Results

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- Experimental Goal: Compare quantification of the 52 carbamates in pure water vs. in waste water
- List of Analyte:

Carbamates Analogs	Thio-Carbamates	Urea Analogs	Complexes & Chelators
Aldicarb Sulfoxide	Butylate	Bramacil	Febam (Fe+3 DMDTC)
Aldicarb Sulfone	Diallate	Chlororoxuran	Metam Na (MMDTC anion)
Aldicarb	DPTC	Diuran	DMDTC (Anion)
Carbaryl	Molinate	Fenuron	OH MMDTC (anion)
Carbofuran	Tillam	Monuron	Ziram (Zn-2 DMDTC)
3OH-Carbofuran	Triallate	Neburon	
Methiocarb	Vernolate	Siduron	
Methomyl		Tebuthiuron	
1-Naphthol			
- System Set-up: Alliance HT 2795/ZQ/PDA990
- Software: QuanOptimize in MassLynx 4.0
- Experimental Conditions:
 - HPLC: Symmetry C8 2.1x150mm, 35µm, 0.3mL/min, 40°C
 - Mobile Phase:
 - A: 10mM ammonium acetate in water at pH 5.0
 - B: 10mM ammonium acetate in acetonitrile
 - Ionization: ESI+
 - Capillary Voltage: 3.5kV
 - Source Temperature: 140°C
 - Desolvation Temperature: 350°C
- Data & Results
 - Standard Spiked into Real Matrices Methomyl

Response type: External Std. Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None
 $r^2 = 0.997705$
0.1 to 100 ppb

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

• Sample Barcode:

• Composition

Ingredient	Compound	Weight (g)
PVC Resin	PV-4410 ☐ Suspension Grade ☒ Emulsion Grade ☐ Bulk Polymerised ☐ Copolymer	2000
Plasticiser	Dinonyl Phthalate (DNP) ▼	700
Stabiliser	Dibasic Lead Phosphite ▼	20
Lubricant	Calcium Stearate ▼	30
Filler	Ground CaCO3 ▼	500
Pigment	DuPont ▼	0

• Processing


1. Weigh all the ingredients and mix them except for PVC resin to a masterbatch.
2. Feed the PVC resin and the masterbatch to a compounding extruder and blend them to a compound with conditions below:
 - Temperature = °C
 - Screw Speed = rpm
3. Cool the blended compound to 60°C and pelletize it to a suitable size.

태그	값
Sample Barcode	19-0001-1-1
PVC Resin	PV-4410
PVC Type	Emulsion Grade
PVC Weight (g)	2000.0
Plasticizer	Dinonyl Phthalate (DNP)
Plasticizer Weight (g)	700.0
Stabilizer	Dibasic Lead Phosphite
Stabilizer Weight	20.0
Lubricant	Calcium Stearate
Lubricant Weight (g)	30.0
Filler	Ground CaCO3
Filler Weight (g)	500.0
Pigment	DuPont
Pigment Weight (g)	0.0
Temperature (°C)	120
Screw Speed (rpm)	300

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태그	ID	값
PVC Weight (g)	PVCWeight	2000.0
Plasticizer Weight (g)	PlasticizerWeight	700.0
Stabilizer Weight	StabilizerWeight	20.0
Lubricant Weight (g)	LubricantWeight	30.0
Filler Weight (g)	FillerWeight	500.0
Pigment Weight (g)	PigmentWeight	0.0
Total Weight (g)	TotalWeight	300

$$\text{TotalWeight} = \text{Sum}(\text{PVCWeight}, \text{PlasticizerWeight}, \text{StabilizerWeight}, \text{LubricantWeight}, \text{FillerWeight}, \text{PigmentWeight})$$

태그	ID	값
Sample Weight (g)	SampleWeight	10.30
Residue Weight (g)	ResidueWeight	9.60
Extractables	Extractables	6.8%

$$\text{Extractables} = (\text{SampleWeight} - \text{ResidueWeight}) / \text{SampleWeight}$$

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



Designation: E 1578 - 06

Standard Guide for

Laboratory Information Management Systems (LIMS)¹

This standard is issued under the fixed designation E 1578; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This guide covers issues commonly encountered at all stages in the life cycle of Laboratory Information Management Systems from inception to retirement. The sub-sections that follow describe details of scope of this document in specific areas.

1.2 *High Level Purpose*—The purpose of this guide includes: (1) help educate new users of Laboratory Information Management Systems (LIMS), (2) provide standard terminology that can be used by LIMS vendors and end users, (3) establish minimum requirements for primary LIMS functions, (4) provide guidance for the specification, evaluation, cost justification, implementation, project management, training, and documentation, and (5) provide an example of a LIMS function checklist.

1.3 *LIMS Definition*—The term Laboratory Information Management Systems (LIMS) describes the class of computer systems designed to manage laboratory information.

1.4 *Laboratory Categories*—The spectrum of laboratories that employ LIMS is wide spread. The following break down provides an overview of the laboratory categories that use LIMS as well as examples of laboratories in each category.

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

The image displays two side-by-side screenshots of the WizLog web application interface, demonstrating its multi-language support. Both screenshots show the 'Write Page' form for a specific note.

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Right Screenshot (English Interface):

- Header:** WizLog®, To Do, Search All, Write Page, View Page, Templates, Binders, 서동건(sdk) logged in at 12:37 PM, Help, Log Out
- Left Sidebar:** Notes, Write Page, Review Page, View Page, Masters, Statistics, User
- Main Content:**
 - Write Page:** Binder (NOTE-002: 고분자 분석), ID (NOTE-002-PAGE-018), No. (18), State (Writing), Title (Quantitative Analysis of Carbamates in Waste Water), References (NOTE-002-PAGE-019: Heavy Metal Analysis of Waste Water [Writing]), File (NOTE-002-PAGE-018 Quantitative Analysis of Carbamates in Waste Water.pdf (187KB)), Created (3/13/2017 12:38:25 PM), Creator (서동건(sdk)), Edited (3/14/2017 10:59:43 AM), Stamped
 - Contents:** Edited (3/14/2017 10:59:43 AM), State (Available), File (T-CHR-001 Automatic LC_MS Quantification.doc (212KB) [Show File]), Add, Edit, Delete, History
 - Experimental Goal:** Compare quantification of the 52 carbamates in pure water vs. in waste water
 - List of Analyte:** Carbamates Analogs, Thio-Carbamates, Urea Analogs, Complexes & Chelators.